10/519,629

GI

```
(FILE 'HOME' ENTERED AT 13:33:01 ON 12 DEC 2005)
```


```
FILE 'REGISTRY' ENTERED AT 13:33:13 ON 12 DEC 2005
L1
               STRUCTURE UPLOADED
L2
             0 S L1
L3
             1 S L1 FULL
     FILE 'CAPLUS' ENTERED AT 13:33:52 ON 12 DEC 2005
             1 S L3
L4
=> d l4 bib abs
     ANSWER 1 OF 1 CAPLUS COPYRIGHT 2005 ACS on STN
L4
     2004:20693 CAPLUS
ΑN
DN
     140:77263
     Preparation of 1,2,4-triazoline-3,5-dione-containing ferrocene compound
ΤI
     and use thereof for determination of vitamin D compound
     Ishigai, Masaki; Murao, Naoaki; Sekiguchi, Nobuo; Takahashi, Tadakatsu
IN
PA
     Chuqai Seiyaku Kabushiki Kaisya, Japan
SO
     PCT Int. Appl., 80 pp.
     CODEN: PIXXD2
DT
     Patent
LA
     Japanese
FAN.CNT 1
                      KIND DATE
                                         APPLICATION NO.
                                                                 DATE
     PATENT NO.
                    71
---- -----
----
                                          -----
     ______
                                                                 -----
                        A1 20040108 WO 2003-JP8166
                                                                20030627
ΡI
     WO 2004002996
        W: AE, AG, AL, AM, AT, AU, AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN,
            CO, CR, CU, CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH,
            GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC, LK, LR,
            LS, LT, LU, LV, MA, MD, MG, MK, MN, MW, MX, MZ, NI, NO, NZ, OM,
            PG, PH, PL, PT, RO, RU, SC, SD, SE, SG, SK, SL, TJ, TM, TN, TR,
            TT, TZ, UA, UG, US, UZ, VC, VN, YU, ZA, ZM, ZW
        RW: GH, GM, KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW, AM, AZ, BY,
            KG, KZ, MD, RU, TJ, TM, AT, BE, BG, CH, CY, CZ, DE, DK, EE, ES,
            FI, FR, GB, GR, HU, IE, IT, LU, MC, NL, PT, RO, SE, SI, SK, TR,
            BF, BJ, CF, CG, CI, CM, GA, GN, GQ, GW, ML, MR, NE, SN, TD, TG
                         A1 20050525 EP 2003-738537
     EP 1533316
                                                                20030627
            AT, BE, CH, DE, DK, ES, FR, GB, GR, IT, LI, LU, NL, SE, MC, PT,
             IE, SI, LT, LV, FI, RO, MK, CY, AL, TR, BG, CZ, EE, HU, SK
PRAI JP 2002-188541
                         Α
                              20020627
     WO 2003-JP8166
                         W
                               20030627
     MARPAT 140:77263
```

## \* STRUCTURE DIAGRAM TOO LARGE FOR DISPLAY - AVAILABLE VIA OFFLINE PRINT \*

Disclosed are a novel ferrocene compound represented by the following formula (I) [Q = a direct bond, alkylene, W1-X-W2; W1 = alkylene, phenylene; W2 = alkylene; X = O, N(Ra)CO, N(Ra)CONH, OCONH, N(Ra)OSO; Ra = lower alkyl; R, R' = H, HO, NO2, cyano, halo, each (un)substituted lower alkyl, lower alkenyl, lower alkynyl, lower alkoxy, lower acyl, CO2H, or CONH2; m = an integer of 1-3; n = an integer of 1-4], a reagent containing the compound, and a high-sensitivity method of determining a vitamin D compound with the reagent. Specifically, the ferrocene compound I is reacted with a VD compound, and the compound comprising these compds. combined with each other, i.e. a Diels-Alder adduct (II) [Q, R, R', m, n = same as above; A1, A3 =(un) substituted lower alkylene, alkenylene, or alkynylene; A2 = a direct bond, CH:CH, C.tplbond.C, O, S, NH; R1 = H, (un)protected OH; R2 = H, HO, halo, (un) substituted lower alkyl, alkenyl, alkynyl, alkoxy, or acyl; R3 = H, a protecting group; R4, R5, R6 = H, H0, N02, cyano, halo, (un) substituted lower alkyl, cycloalkyl, lower alkenyl, lower alkynyl, lower alkoxy, lower acyl, CO2H, CONH2, or NH2; R7, R8 = H, OH; or R7 and R8 together form a double bond] is subjected to LC/ESI-MS/MS. Thus, the VD compound can be determined with higher sensitivity than in conventional

techniques. The ferrocene compound is extremely useful as an agent for derivative formation when a VD compound is determined by LC/ESI-MS/MS. The compound obtained, which comprises the ferrocene compound and VD compound which have been combined with each other, is useful as, e.g., a labeled compound in the determination of a VD compound by LC/ESI-MS/MS. Alfacalcidol in rat serum was determined by adduct formation with 4-(ferrocenylmethyl)-1,2,4-triazoline-3,5-dione at the lower detection limit of 0.08 ng/mL which was 125-times more sensitive than that of the direct method (10 ng/mL).

RE.CNT 2 THERE ARE 2 CITED REFERENCES AVAILABLE FOR THIS RECORD ALL CITATIONS AVAILABLE IN THE RE FORMAT



10/519,629

Ref #	Hits	Search Query	DBs	Default Operator	Plurals	Time Stamp
L1	99	(548/103).CCLS.	US-PGPUB; USPAT; EPO; JPO	OR	OFF	2005/12/12 14:05
L2	177	(556/144).CCLS.	US-PGPUB; USPAT; EPO; JPO	OR	OFF	2005/12/12 14:12
L3	872	(514/167).CCLS.	US-PGPUB; USPAT; EPO; JPO	OR	OFF	2005/12/12 14:28
L4	340	(514/502).CCLS.	US-PGPUB; USPAT; EPO; JPO	OR	OFF	2005/12/12 14:28